


Leilox

1.5/2.0 MicroLocking System



Special Thanks to Dr. Hugo Schmökel

R **RITA**
LEIBINGER
MEDICAL



LeiLOX 1.5 / 2.0

Titanium Locking Fracture System for cats and small dogs

Accidents happen. And you, as a surgeon, are often required to fix complex fractures.

Our 1.5/2.0 locking plate system features 32 different plates, making it the most flexible titanium locking system for small animals. The system can handle both simple and complex fractures.



Titanium: Best Biocompatibility

Unique Titanium Alloy makes our implants very strong, resistant and lightweight, which is ideal for small animals. Titanium offers advantages in biocompatibility, extreme temperatures and weight.



All-in-One Set

LeiLOX 1.5/2.0 is the "swiss army knife" for micro fractures: it contains not only (cuttable) standard plates, but also unique shapes for complex anatomy, in a compact instrument set.



Various Applications

With more than 32 different plates, LeiLOX 1.5/2.0 is the perfect solution to treat a wide variety of fracture cases: Y-Plates, L-Plates, Bridge-Plates, Standard-Plates and Cuttable 20-hole Plates offer plenty of possibilities: from maxillofacial and pelvic fractures, to the fixation of long tubular bones.

Stardrive Screws

All of our locking screwheads are stardrive, which allows for better tightening torque than hexagonal screwheads.

Additionally, you do not need a separate holding sleeve if you use a screwdriver shaft from Rita Leibinger. The screwdriver shaft is self-holding and keeps the screws firmly on the driver.



BEST QUALITY
MADE IN GERMANY

Osteosynthesis Locking Set 1.5/2.0

Osteosynthese Locking Set 1.5/2.0 Titan


Contains

- 1 locking plate of each style
- 5 screws of each Length
- 2 star drive screwdriver shafts
- 1 screwdriver handle
- 4 locking drill guides of each diameter
- 1 depth gauge
- 2 drills (AO)
- 1 pair of bending irons.
- Tray inclusive.

142-1520-00

Tray without contents.

142-1520-10



Sterilization Tray

for Trays, 310x190x130mm
blue (image)

150-5401-30

green

150-5402-30



1520.leibinger.vet




also Online available
veterinary.shop

LeiLOX Starter Set 1.5

Contains:

- 5 screws of each Length
- Bone Plate Locking, 20 hole, cuttable
- Bone Plate Locking, 4 Hole, 14 mm Bridge
- Bone Plate Locking, 6 Hole, 21 mm Bridge
- Bone Plate Locking, 6 Hole, 35 mm Bridge
- Locking-L-Plate left, Locking-L-Plate right
- Locking-Y-Plate
- Locking drill guide
- Screwdriver shaft Star Drive
- Drill AO
- Bone Plate Cutter

142-1501-00




LeiLOX Starter Set 2.0

Contains:

- 5 screws of each Length
- Bone Plate Locking, 20 hole, cuttable
- Bone Plate Locking, 4 Hole, 18 mm Bridge
- Bone Plate Locking, 6 Hole, 27 mm Bridge
- Bone Plate Locking, 6 Hole, 36 mm Bridge
- Locking-L-Plate left, Locking-L-Plate right
- Locking-Y-Plate
- Locking drill guide
- Screwdriver shaft Star Drive
- Drill AO
- Bone Plate Cutter

142-2001-00



Locking Plates 1.5

1.5 Locking Plate Titanium

5mm broad, 1.5mm thick



Product Code	Description
142-0215-04	4-Hole, 25.5mm
142-0215-05	5-Hole, 32.5mm
142-0215-06	6-Hole, 39.5mm
142-0215-07	7-Hole, 46.5mm
142-0215-08	8-Hole, 53.5mm
142-0215-09	9-Hole, 60.5mm
142-0215-10	10-Hole, 67.5mm
142-0215-20	20-Hole, cuttable

1.5 Locking Plate with Bridge, Titanium

5mm broad, 1.5mm thick



Product Code	Description
142-0215-614	14mm Bridge, 4-Hole, 32.5mm
142-0215-621	21mm Bridge, 6-Hole, 53.5mm
142-0215-635	35mm Bridge, 6-Hole, 67.5mm
142-0215-658	58mm Bridge, 6-Hole, 95mm

1.5 Locking-L-Plate Titanium

5mm broad, 1.5mm thick, 25.5 x 11.5mm Length

Product Code	Description
142-0215-504L	left
142-0215-504R	right

1.5 LeiLOX Locking Screw Titan

for LeiLox Locking System, Star Drive T6 self-holding (T6 Shaft from Rita Leibinger needed) self-tapping with three flute cutting edge



Length (mm)	Product Code
5	245-415-05
6	245-415-06
7	245-415-07
8	245-415-08
9	245-415-09
10	245-415-10
11	245-415-11
12	245-415-12
14	245-415-14
16	245-415-16
18	245-415-18
20	245-415-20
22	245-415-22
24	245-415-24
26	245-415-26

1.5 Locking-Y-Plate Titanium

5mm & 11mm broad, 1.5mm thick

Product Code	Description
142-0215-805	5-Hole, 33mm long
142-0215-809	10-Hole, 61mm long



2.0 Locking Plate Titanium

6mm broad, 1.5mm thick



Product Code	Description
142-0220-04	4-Hole, 33mm
142-0220-05	5-Hole, 42mm
142-0220-06	6-Hole, 51mm
142-0220-07	7-Hole, 60mm
142-0220-08	8-Hole, 69mm
142-0220-09	9-Hole, 78mm
142-0220-10	10-Hole, 87mm
142-0220-20	20-Hole, cuttable

2.0 LeiLOX Locking Screw Titan

for LeiLox Locking System, Star Drive T8 self-holding (T8 Shaft from Rita Leibinger needed) self-tapping with three flute cutting edge



Length (mm)	Product Code
5	245-420-05
6	245-420-06
7	245-420-07
8	245-420-08
9	245-420-09
10	245-420-10
11	245-420-11
12	245-420-12
14	245-420-14
16	245-420-16
18	245-420-18
20	245-420-20
22	245-420-22
24	245-420-24
26	245-420-26
28	245-420-28
30	245-420-30



2.0 Locking Plate with Bridge, Titanium

6mm broad, 1.5mm thick



Product Code	Beschreibung
142-0220-618	18mm Bridge, 4-Hole, 42mm
142-0220-627	27mm Bridge, 6-Hole, 69mm
142-0220-636	36mm Bridge, 6-Hole, 78mm
142-0220-668	68mm Bridge, 6-Hole, 110mm

2.0 Locking-L-Plate Titanium

6mm broad, 1.5mm thick, 33 x 15mm Length

Product Code	Description
142-0220-504L	left
142-0220-504R	right



2.0 Locking-Y-Plate Titanium

6mm & 13.5mm broad, 1.5mm thick

Product Code	Description
142-0220-805	5-Hole, 42mm long
142-0220-809	10-Hole, 78mm long

